



OEM Peristaltic Pump Without Drive

T-S500&WX10-14-H



- Use WX10-14-H pump head and step motor.
- Ideal OEM peristaltic pump for equipment and instruments.
- Soleplate installation is standard, panel installation is optional.
- Meet RoHS requirement.



Pump Specification

Speed	≤100rpm
Max. Flow Rate	40mL/min
Tubing	ID≤3.17mm, wall thickness 0.8mm-1.0mm
Dimension (L×W×H)	106×61×81 (mm)
Relative humidity	<80%
Operating temperature	0 to 40°C

Motor Specifications


Step angle	1.8°
Phase	2
Phase voltage	6.56V
Phase current	0.8A
Resistance	8.2Ω±10%
Inductance	13.9mH±20%

Motor Port


Line Color	Definition
Brown	A
Orange	\bar{A}
Red	B
Yellow	\bar{B}

OEM Peristaltic Pump Without Drive

T-S502&JY10-14-H



- Use JY10-14 pump head and step motor.
- Ideal OEM peristaltic pump for equipment and instruments.
- Panel installation is standard, Soleplate installation is optional.
- Meet RoHS requirement.



Pump Specification

Speed	≤60rpm
Max. flow rate	24mL/min
Tubing	ID≤3.17mm Wall thickness 0.8-1.0(mm)
Dimension(L×W×H)	66×82×100(mm)
Relative humidity	<80%
Operating temperature	0 to 40°C

Motor Specifications

Step angle	1.8°
Phase	2
Phase voltage	8.4V
Phase current	0.6A
Resistance	14Ω±10%
Inductance	23.5mH±20%

Motor Port

Line Color	Definition
Red	A
Green	\bar{A}
Yellow	B
Blue	\bar{B}